

SuperEar®

**Assistive Listening Device
ADA Model SE5200A**

Manufacturer's Suggested Cleaning Recommendations



The SE5200A is packaged with one pair of foam ear pads in place on the headphones, and three additional pairs of foam pads may be included beneath the packing tray. Additional quantities of foam ear pads may be purchased from the manufacturer at \$0.50 per pair.

It is recommended that, after each use by a patient, the foam ear pads be removed and discarded. If desired, the earphones may be cleaned with a light application of an alcohol-based disinfectant before replacing with new foam ear pads.

SonicTechnology®

120 Richardson St., Grass Valley, CA 95945
(530) 272-4607 ■ (800) 247-5548 ■ FAX (530) 272-4257
www.sonictechnology.com